

Tecnis Multifocal Lenses for Cataract Surgery

Offers you the option to choose the right lens for your lifestyle. Whether you want more comfort while reading a book, more confidence viewing your smartphone, or spend much time on a laptop or desktop, the family of TECNIS Multifocal IOLs provides just that. With the TECNIS Multifocal, in addition to regaining their distance vision, patients with age-related vision loss can look forward to regaining their near vision and reducing their dependence on glasses based on the three options for near vision below.

ZMB00 1-Piece +4.0 Near Add

The ZMB00 is optimized for patients favoring up close near-vision activities such as reading and knitting. Enhanced near-vision performance for this IOL is at ~33 cm (13 inches) theoretical reading distance.



ZLB00 1-Piece +3.25 Low Add

The ZLB00 is optimized for most patients' average reading distance, favoring activities such as newspaper and book reading, as well as use of a phone and tablet. Enhanced near-vision performance for this IOL is at ~42 cm (16.5 inches) theoretical reading distance.



ZKB00 1-Piece +2.75 Low Add

The ZKB00 is optimized for patients wanting a balance between reading and intermediate-vision activities such as lap reading, computer vision, and other intermediate vision tasks such as reading the dashboard on their car and grocery shopping. Enhanced near-vision performance for this IOL is at ~50 cm (19.5 inches) theoretical reading distance.

